

RECOMMENDED AND MAXIMUM ROTATION SPEEDS

Recommended speed: 7.000 – 10.000 min
for Coarse/Medium/Fine

Maximum speed: 20.000 rpm

KATANA™ Zirconia TWIST DIA

Zirconia finishing and polishing system

The specifications and appearance of the product are subject to change without notice. Printed color can be slightly different from actual color. "KATANA" and "CERABIEN" are trademarks of Noritake Co., Ltd. "PANAVIA" and "CLEARFIL" are trademarks of Kuraray Co., Ltd. "CEREC" and "SpeedFire" are trademarks of Dentsply Sirona Inc.

Kuraray Europe GmbH, Philipp-Reis-Str. 4 65795 Hattersheim am Main, Deutschland

Phone +49 (0) 69-305 85 980 E-Mail dental@kuraray.com



Website www.kuraraynoritake.eu



POLISHING ZIRCONIA WITH KATANA™ Zirconia TWIST DIA

Fine gloss in 3 steps

Flexible diamond impregnated spirals for polishing zirconia.



WONDERFUL HIGH-SHINE RESULTS ON ZIRCONIA

- Flexible polishing twist for refining all types of tooth surfaces.
- Highly suitable for occlusions
- Modelling of zirconia restoration is maintained
- Sterilizable and reusable
- Economic efficiency due to high durability and reusability.




Occlusal



Buccal



Incisal



Labial



Lingual



Approximal