To meet your need for perfection

... all things cotton

.....

#-

h

гоеко

- Reliable high manufacturing standard qualityfocussed and fully automated
- Pharmaceutical quality raw materials tested for purity, foreign fibres, absorbency and surfactant residues
- Versatile and well thought out adapted to the needs of dentists
 - Cotton products Made in Germany

#COLTENE

roeko – one name – long-standing experience – extensive technical know-how – utmost attention to quality

The quality of cotton products definitely pays dividends if one does the sums correctly. The reliable high manufacturing standard of ROEKO cotton products contributes to lower consumption,

trouble-free working and hence to truly perfect results.



LUNA GOLD

See and feel the difference

LUNA GOLD is encased in skin-friendly mesh and can be removed without prior wetting – this prevents mucosal reaction and saves time. The gold stripe clearly differentiates LUNA GOLD from conventional dental rolls.

For regular use we recommend

LUNA – long-lasting absorbency and exceptional dimensional stability

COTTON-ROLL – particularly soft with rapid absorption.

ROEKO dental rolls are bleached chlorine-free.

ROLL-O-MAT

The system is both simple and ingenious: slight pressure on the lid is sufficient to dispense a single cotton roll and keep it positioned reliably on the dispenser tip. The ROLL-O-MAT is available as plastic dispenser or in a stainless steel version.





LUNA FOR KIDS

Specially designed for the very small patients LUNA FOR KIDS dental rolls are particularly

short and yet dimensionally stable like their larger Luna counterparts. This leaves plenty of room to work, even in the smallest of mouths. ROEKO dental rolls are bleached chlorine-free.

LUNAMAT

The LUNAMAT is equipped with a drawer which dispenses 2-3 cotton rolls when pulled out.



PAROTISROLL

Multipurpose cotton roll with spring action

The spring rod inside the cotton roll provides a spring action and makes the PAROTISROLL self-fixating. When used buccally, lingually or labially, it keeps the mouth open, holds down the tongue and retracts the lips. When positioned appropriately, it also prevents the flow of saliva through pressure on the parotid gland.

PAROTISROLL ANATOMIC

By Dr. Hannelore Oberschachtsiek PAROTISROLL anatomic is the Parotisroll with flattened ends. The flattened ends facilitate placement of the cotton roll and reduce the foreign body feeling in the mouth.

PAROTIMAT

For convenient individual dispensing and hygienic storage of 50 PAROTISROLLS.





ROEKO dental roll dispensers keep cotton rolls within reach during treatment while at the same time protecting them hygienically from spray and dust.



COTTON PELLETS

Uniform in shape and size

The fully automated production of the ROEKO cotton pellets guarantees their uniform shape and size and thus, reliability. They are the ideal aid for drying cavities or applying medications.

#000	#00	#0	#1	#2
0	\bigcirc	\bigcirc	\bigcirc	\bigcirc

Ø3mm Ø4mm Ø4,8mm Ø6,3mm Ø7,6mm

SOLOMAT

The elegant dispenser

In this stainless steel dispenser, the cotton pellets are spring-loaded to the dispensing surface. The scratch- and abrasion-resistant surface prevents bacteria from adhering and is easy to clean and disinfect.



SOLOMAT-N

For individual dispensing of cotton pellets or foam pellets

Autoclavable plastic dispenser with refill system. The cotton pellets are moved manually to the dispensing surface.



COMPRECAP

COMPRECAP compression caps are a simple and intelligent aid for preparing impression taking. They stop bleeding – by compression and dry tooth stumps.

COMPRECAP ANATOMIC

COMPRECAP ANATOMIC are ideally adapted to the anatomy of the dental arch through semi-circular cut-outs.

SOLOMAT-CLICK

The smart dispensing system SOLOMAT-CLICK consists of a dispensing grille, cover and refill. These 3 parts are clicked together when repurchasing a refill. The stainless steel base gives the dispenser the necessary stability.



The ROEKO cotton pellet dispensers keep cotton pellets within reach during treatment while at the same time protecting them hygienically from spray and dust. A firm base and well-designed lid allow dispensing cotton pellets with only one hand. The rounded rods of the dispensing area prevent the pellets from being spueezed on removal as opposed to a dispensing area with a grid structure.

Order information

	Contents	REF
LUNA	300 g	
		100 001
		100 002
		100 003
	# 4	100 004
	500 pcs.*	
		105 001
		105 002
	# 3	105 003
	1000 pcs.*	
		104 001
		104 002
		104 003
LUNA FOR KIDS	200 pcs.	
	# 2	100 032
LUNA GOLD	300 g	
		109 121
		109 122
	# 3	109 123
	1000 pcs.*	
		109 101
		109 102
	# 3	109 103
COTTON-ROLL	1000 pcs.*	
		104 011
		104 012
	# 3	104 013
ROLL-O-MAT	1 pc.	
plastic, white		120 004
$134 \times 100 \times 48 \text{ mm}$		
chrome nickel steel		120 005
134 × 100 × 43,5 mm (can be		
used with cotton roll sizes 1–3), standard length (36 – 40 mm). Not		
suitable for LUNA FOR KIDS.		
LUNAMAT	1 рс.	
plastic, white		120 000
$52 \times 106 \times 100$ mm (can be used with cotton roll sizes 1– 4), stan-		
dard length $(36 - 40 \text{ mm})$ and with		
LUNA FOR KIDS.		

SOLOMAT- CLICK ² plastic 43 × 50 × 41 mm		# 000 # 00 # 0 # 1	1 pc.	224 000 224 100 224 110 224 111
REFILLS		# 000 # 00 # 0 # 1	3 pcs.	214 000 214 100 214 110 214 111
SOLOMAT-N ² plastic 75 × 43 × 52 mm		# 000 # 00 # 1 # 2	1 pc.	223 000 223 100 223 110 223 111 223 111 223 112
REFILLS		# 000 # 00 # 1 # 2	6 pcs.	213 000 213 100 213 110 213 111 213 111 213 112
COMPRECAP			120 pcs. 600 pcs. 120 pcs. 600 pcs. 60 pcs. 400 pcs. 120 pcs.	530 008 531 008 530 010 531 010 530 014 531 014 530 000
COMPRECAP ANATOMIC		# 1 # 3 # 5 Set (#1/3/5)	120 pcs. 120 pcs. 60 pcs. 120 pcs.	531 201 531 203 531 205 531 299
Dispenser is supplied with	nout cotton pe	ellets		

Contents

¹⁾ Dispenser is supplied without cotton pellets

²⁾ Dispenser is supplied with cotton pellets

COTTON PELLETS

SOLOMAT¹

Coltène/Whaledent	Coltène/Whaledent Inc.

235 Ascot Parkway Cuyahoga Falls, Ohio 44223 / USA T +1 330 916 8800 F +1 330 916 7077 info.us@coltene.com



Coltène/Whaledent Ltd.

The President Suite Kendal House, Victoria Way Burgess Hill, West Sussex RH15 9NF / UK T +44 1444 235486 F +44 1444 870640 info.uk@coltene.com

© 2016 Coltène/Whaledent AG - www.coltene.com

* Only for export

GmbH + Co. KG Raiffeisenstraße 30 89129 Langenau / Germany T +49 7345 805 0 F +49 7345 805 201

info.de@coltene.com