



HAHNENKRATT.COM

PROVISA

THE ECONOMICAL TEMPORARY SOLUTION



Aluminium



Tin

Individual adaption

by using a shell and modulation of the fissures by the occlusion

Free of nickel

For patients with metallic dental fillings (e.g. Amalgam, Gold) in their dentures, we recommend to use PROTECT Crowns REF 480X... with passivated, gold coloured oxide layer.

Economical and efficient

- time-saving, because there is no need to make an impression
- cost-saving, because of low material costs

Optimal fixation of the occlusion

No harmful influence to the vitality of the pulpa

While using temporary composites, high temperature can occur which can be harmful for the pulpa.*)

Tested concerning biocompatibility according to ISO 10993

*) Literature: Scientific report in Quintessenz 51,4,349-356(2000)



PROVISA Shells	REF		
	Contents	Aluminium Shells	Tin Shells
Plastic assortment	100	372	378
Cardboard assortment	100	371	
Refill, Please add size	100	370X...	
Refill, Please add size	10	374X...	382X...